

Comfort-4 Eco

THERAPY CHAIR

HIGHEST COMFORT, ATTRACTIVE DESIGN IN ADDITION TO FUNCTIONALITY AND SAFETY

- The Comfort-4 Therapy Chairs are designed to assist and facilitate the work of physicians and medical staff
- They increase the patient's comfort considerably and therefore help to relieve tension during treatment
- Several optional accessories for optimal choice



PROFESSIONAL SEATING SOLUTIONS BY

Digiterm[®]
medical equipment

* Chair in the photo involves some optional elements

Comfort-4 Eco

ART.-NR. 723C4E00

Optional accessories	Art.-Nr.
EPN kit incl. multicontact pin with internal connecting cable HIGHLY RECOMMENDED FOR DIALYSIS	713E004
Trendelenburg foot-switches on both sides of underframe	713E005
Rechargeable battery back-up operation in case power failure (24V DC)	713E006
Reading lamp with LED light	713E008
USB charging unit for smart devices (5V, 1A)	713E009
Connecting cable (3m) for potential equalization (EPN)	713E012
Comfort-T2 Multifunctional table	713M009
Flexible height adjustable armrest for universal medical applications (one pair)	713M010
Integrated gas spring in the arm support for smoother and safer operation (one pair)	713M011
Push bar	713M012
Universal holder for tablets and smartphones	713M013
Handle on both sides of the seat	713M014
Removable footrest for standalone exercise applications	713M031
Paper roll holder	713M033
Motorized footrest kit (surcharge of upgrade from manual footrest)	713M036
Infusion bottle holder (IV pole) with height adjustment	713M038
Central operating brake with 4x diam. 10 cm twin wheel castors	713M040
Central operating total and directional lock with 4x diameter: 12.5 cm twin wheel castors	713M041
Individual logo on the upholstery	713U011
Headrest pillow square	713U012
DuoShape upholstery in individual colour combination (40 diff. colours) TUNDRA	713U026
Transparent protecting plastic cover for the footrest	713U027
Transparent protecting plastic cover for the legrest adjustable between 0° and -30° below horizontal	713U028
DuoShape upholstery with 110 mm thickness in standard Chrom-Perle colour combination	713U031
DuoShape upholstery with 110 mm thickness in individual colour combination (40 diff. colours) TUNDRA	713U032
DuoShape upholstery with memory foam in standard Chrom-Perle colour combination	713U033
DuoShape upholstery with memory foam in individual colour combination (40 diff. colours)	713U034
DuoShape upholstery with flame retardant artificial leather combination (15 diff. colours) PALMA	713U035

General technical data

Seat height adjustment	55–77 cm	Dimensions of armrest	L 57 × W 15 × H 6,5 cm
Backrest adjusting range	80° (-13° under horizontal in Trendelenburg position)	Size of castors	4 × Ø 10 cm swivel castors with individual total lock
Seat tilt adjustment	20°	Weight of chair	80 kg
Legrest adjustment	30°	Maximal load	200 kg weight capacity
Footrest adjustment	25 cm (manually)	Standard Power Supply	100–240 V AC, 50/60 Hz
Dimensions in seat position	L 145 × W 92.5 or L 145 × W 80.5 cm*	Technical data of motors	24V DC, IP X6
Size of upholstery	L 210 × W 60 × H 9.5 cm	Applied directive and standards	93/42/EEC and harmonized standards (EN60601-1; EN60601-2-52; etc.)

* Depend on armrest and footrest setting



Digiterm Ltd.
Ipar út 7, H-3200 Gyöngyös
Tel.: +36 37 311 990
Fax: +36 37 311 893
E-mail: sales@digiterm.hu
www.digiterm.hu

