

Comfort-4 Battery THERAPY CHAIR

HIGHEST COMFORT, ATTRACTIVE DESIGN IN ADDITION TO FUNCTIONALITY AND SAFETY

- The Comfort-4 Therapy Chairs are designed to assist and facilitate the work of physicians and medical staff
- They increase the patient's comfort considerably and therefore help to relieve tension during treatment
- Several optional accessories for optimal choice



PROFESSIONAL SEATING SOLUTIONS BY

Digiterm[®]
medical equipment

* Chair in the photo involves some optional elements

Comfort-4 Battery

ART.-NR. 723C4B00

Optional accessories	Art.-Nr.
EPN kit incl. multicontact pin with internal connecting cable HIGHLY RECOMMENDED FOR DIALYSIS	713E004
Trendelenburg foot-switches on both sides of underframe	713E005
Connecting cable (3m) for potential equalization (EPN)	713E012
Reading lamp with flexible arm	713E013
Comfort-T2 Multifunctional table	713M009
Flexible height adjustable armrest for universal medical applications (one pair)	713M010
Integrated gas spring in the arm support for smoother and safer operation (one pair)	713M011
Push bar	713M012
Universal holder for tablets and smartphones	713M013
Handle on both sides of the seat	713M014
Removable footrest for standalone exercise applications	713M031
Paper roll holder	713M033
Infusion bottle holder (IV pole) with height adjustment	713M038
Central operating brake with 4x diam. 10 cm twin wheel castors	713M040
Central operating total and directional lock with 4x diameter: 12.5 cm twin wheel castors	713M041
Stainless steel armrest protecting frame (one pair)(together with armrest longitudinal adjustment)	713M045
Individual logo on the upholstery	713U011
Headrest pillow square	713U012
DuoShape upholstery in individual colour combination (40 diff. colours) TUNDRA	713U026
Transparent protecting plastic cover for the footrest	713U027
Transparent protecting plastic cover for the legrest adjustable between 0°and -30° below horizontal	713U028
DuoShape upholstery with 110 mm thickness in standard Chrom-Perle colour combination	713U031
DuoShape upholstery with 110 mm thickness in individual colour combination (40 diff. colours) TUNDRA	713U032
DuoShape upholstery with memory foam in standard Chrom-Perle colour combination	713U033
DuoShape upholstery with memory foam in individual colour combination (40 diff. colours) TUNDRA	713U034
DuoShape upholstery with flame retardant artificial leather combination (15 diff. colours) PALMA	713U035

General technical data			
Seat height adjustment	55–77 cm	Size of castors	4x Ø 10 cm swivel castors with individual total lock
Backrest adjusting range	80° (-13° under horizontal in Trendelenburg)	Net weight of chair	80 kg
Seat tilt adjustment	20°	Maximal load	200 kg weight capacity
Legrest adjustment	30°	Technical data of motors	24V DC, IPX6
Footrest adjustment	25 cm (manually)	Standard Power Supply	230 V, 50/60 Hz (100 V, 120 V optional)
Size of upholstery	L 210 × W 60 × H 9.5 cm	Rechargeable battery back-up	24 V/1,2 Ah
Armrest size	L 57 × W 15 × H 6.5 cm	Applied directive and standards	93/42/ EEC and harmonized standards (EN60601-1; EN60601-2-52; etc.)



Digiterm Ltd.
Ipar út 7., H-3200 Gyöngyös
Tel.: +36 37 311 990
Fax: +36 37 311 893
E-mail: sales@digiterm.hu
www.digiterm.hu

